

This Page Is Inserted by IFW Operations  
and is not a part of the Official Record

## **BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images may include (but are not limited to):

- BLACK BORDERS
- TEXT CUT OFF AT TOP, BOTTOM OR SIDES
- FADED TEXT
- ILLEGIBLE TEXT
- SKEWED/SLANTED IMAGES
- COLORED PHOTOS
- BLACK OR VERY BLACK AND WHITE DARK PHOTOS
- GRAY SCALE DOCUMENTS

**IMAGES ARE BEST AVAILABLE COPY.**

**As rescanning documents *will not* correct images,  
please do not report the images to the  
Image Problem Mailbox.**

①



Europäisches Patentamt  
European Patent Office  
Office européen des brevets

⑪ Publication number:

**0 232 706**  
**A2**

⑫

# EUROPEAN PATENT APPLICATION

⑲ Application number: 87100099.8

⑳ Date of filing: 07.01.87

⑤① Int. Cl.: **G 01 N 33/80, G 01 N 33/52,**  
**G 01 N 33/574**  
**// G01N33/577**

③① Priority: 10.01.86 US 817874

④③ Date of publication of application: 19.08.87  
Bulletin 87/34

⑥④ Designated Contracting States: DE FR GB

⑦① Applicant: Sloan-Kettering Institute For Cancer  
Research, 1275 York Avenue, New York New  
York 10021 (US)

⑦② Inventor: Cordon-Cardo, Carlos, 430 East 67th Street,  
New York New York 10021 (US)  
Inventor: Lloyd, Kenneth O., 4525 Henry Hudson  
Parkway West, Bronx New York 10471 (US)  
Inventor: Finstad, Connie L., 504 East 81st Street, New  
York New York 10028 (US)  
Inventor: McGroarty, M.E., 16 Franklin Place, Great  
Neck New York 10021 (US)  
Inventor: Bander, Neil H., 1385 York Avenue, New York  
New York 10021 (US)  
Inventor: Old, Lloyd J., 600 West End Avenue, New York  
New York 10024 (US)  
Inventor: Melamed, Myron R., 35 Braeside Lane, Dobbs  
Ferry New York 10522 (US)

⑦④ Representative: Patentanwälte Schulze Horn und  
Hoffmeister, Goldstrasse 36, D-4400 Münster (DE)

⑤④ Blood group antigen panel.

⑤⑦ Monoclonal antibodies specific for human blood group  
antigens are useful in cancer diagnosis. It has been found that  
cancer patients often either express antigens of a blood type  
different from their normal blood type or cease to produce nor-  
mal blood group antigens in cancerous tissues. By examining  
tissue samples with the monoclonal antibodies described here-  
in, it becomes possible to diagnose cancer.

**EP 0 232 706 A2**

1

FIELD OF THE INVENTION

5           This invention relates to the use of monoclonal  
antibodies in determining the presence of particular  
antigens. Applications of this invention include, but are  
not limited to, diagnosis of disease, including cancer cell  
10       typing or classification and identification of precancerous  
lesions.

PRIOR ART

15           Blood group antigens are carbohydrate determinants  
which are typically found on erythrocytes, certain  
epithelial tissues, and in body secretions. They are formed  
by the sequential addition of saccharides to carbohydrate  
20       side chains of lipids and proteins. Hakomori, Seminars in  
Hematology 18:39 (1981). Genes control synthesis of these  
structures, as well as their expression in secretions, and  
on cell types other than erythrocytes (red blood cells).

25           The "A", "B" and "H" blood group antigens are  
known, at least indirectly, as identifying blood "type".  
Presence of "H" antigen only is characteristic of "Type O"  
30       blood, whereas presence of antigen "A" and "B" in the same  
sample is characteristic of type "AB". Presence of "A" or  
"B" antigen is characteristic of, respectively, type A or  
type B blood. Lewis antigens, i.e., Le<sup>a</sup> and Le<sup>b</sup>, are  
35       typically found in plasma, secretions, and secretory  
epithelia.

1

These antigens are characteristic of other  
5 conditions, and determination of their presence is useful in  
areas other than blood typing. For example, Emmott, et al.,  
J. Urol. 121:37 (1979) have found that the antigens of the  
ABH system which are usually present in normal urinary  
10 bladder tissue, are absent in urinary bladder tumors.  
Additional studies have shown that loss of these antigens is  
an early event in malignant transformation. Liss, et al.,  
Am. J. Clin. Pathol. 68:372(1977); (larynx carcinoma);  
15 Weinstein, et al., Cancer 43:661 (1979) (urinary bladder  
carcinoma). In patients with epithelial cancers, especially  
colon carcinomas, elevated levels of Le<sup>a</sup> and Le<sup>b</sup> antigens  
have been found. Koprowski, et al., Science 212:53 (1981).  
20 Additionally, the presence of normally incompatible blood  
groups in the same patient has been described in some cancer  
patients. Hatton, et al., Biochem. Biophys. Acta 666:361  
(1981).

25

Determination of a change in the amount of blood  
group antigen, a sudden appearance, or disappearance, is  
indicative of a pathological state. Hatton, et al., supra,  
30 for example, found A antigen in type O cancer patients.  
Emmol, Lin, and Koprowski, supra have all shown that  
disappearance or appearance of antigens is typical of  
cancer. Hence, it is desirable to have a method for  
35 determining the presence of blood group antigens.

1

5               Several monoclonal antibodies, which are known to  
be specific for blood group antigens have been used to form  
a "panel" for determining particular antigens. H29-36  
monoclonal antibody, for example, determines the presence of  
10 all varieties of A antigen. See, e.g., Sakamoto, et al.,  
(unpublished manuscript), copending U.S. Patent Application  
Serial No. 474,415, "Monoclonal antibody S8 is known to  
detect B-antigen," Ueda, et al., PNAS 78:5122 (1981).  
15 Additionally, monoclonal antibodies T-174, T-218, P-12, and  
F-3 are specific for Le<sup>a</sup>, Le<sup>b</sup>, X, and Y antigens,  
respectively. Antibody K-21 detects precursor type antigen.  
Rettig, et al., Cancer Res. 45:815 (1985), Lloyd, et al.,  
20 Immunogenetics 17:537 (1983).

              These monoclonal antibodies all of which are  
described more fully, infra, are used as part of a panel to  
25 determine blood group antigens.

              Further details on the panel, its uses, and  
further embodiments, are presented in the description which  
30 follows.

#### BRIEF DESCRIPTION OF THE FIGURES

              Figure 1 shows pictorially the structures and  
35 origins of Le<sup>a</sup>, Le<sup>b</sup>, H-1, H-2, X, Y, A and B antigens.

1

5

Figures 2 and 3 illustrate immunohistological staining patterns of monoclonal antibodies described herein, when applied to normal human adult kidney and urothelium.

10

DETAILED DESCRIPTION OF PREFERRED EMBODIMENTS

Antibodies

15

The hybridoma cell lines which produce the monoclonal antibodies of this invention have been deposited with the American Type Culture Collection, 12301 Parklawn Drive, Rockville Maryland 20852 and bear the following accession numbers:

20

	<u>Hybridoma</u>	<u>ATCC #</u>
	H 29-36	HB 8248
	S 8	
25	T 174	HB 8242
	T 218	HB 8249
	P 12	HB 8551
	F 3	HB 8217
30	K 21	HB 8549

35

Information on derivation of these hybridomas may be found in copending U.S. Application Serial No. 474,415 (H 29-36, T-174, T-218), Serial No. 297,814 (S-8); Serial No.

1

5 604,080 (P-12) and K-21) and Serial No. 470,815 (F-3). In  
addition, the hybridomas are described in Ueda, et al., PNAS  
78:5122 (1981) (S8); Rettig, et al., Cancer Res. 45:815  
10 (1985) (P-12, and K-21) and Lloyd, et al., Immunogenetics  
17:537 (1983) (F-3). The disclosures of all of these are  
incorporated by reference herein.

In summary, the hybridomas are prepared following  
15 the Kohler-Millstein method well known to the art, using, as  
immunizing cell lines, the materials set forth in Table 1.

20

25

30

35

0232706

1107-01-87

SK 375

Table 1. Derivation and specificity of mouse monoclonal antibodies identifying Blood Group Antigens.

Antibody (Ig subclass)	Immunizing Cell Line	Blood Group Specificity.
K21 ( $\mu$ )	Tera-1 Teratocarcinoma	Precursor (Type 1 Chain)
T174 ( $\gamma$ 1)	SK-CO-10 Colon Cancer	Le <sup>a</sup> (Type 1 Chain)
T218 ( $\mu$ )	SK-CO-10 Colon Cancer	Le <sup>b</sup> (Type 1 Chain)
P12 ( $\mu$ )	Fresh Human Placenta	X (Type 2 Chain)
F3 ( $\mu$ )	SK-LU-3 Lung Cancer	Y (Type 2 Chain)
T36 ( $\gamma$ 3)	HT29 Colon Cancer	A (Type 1 and 2 Chains)
S8 ( $\mu$ )	SK-RC-7 Renal Cancer	B (Type 2 Chain)

Note: Purified agglutinin I from Ulex Europeans (Vector Laboratories, Burlingame, CA) served to identify the H-antigen.



1

5 Because expression of ABO(H) antigens on the urogenital system has been discussed in the literature, e.g., Coon, et al., Am. J. Clin. Path. 76:163 (1981); Szulman, J. Exp. Med. 111:785 (1960), and recent studies

10 have shown localization of Le<sup>a</sup> and Le<sup>b</sup> antigens in normal adult urothelium; Juhl, J. Histochem. Cytochem. 33:309 (1985), the urinary tract was used as an exemplary system. The following experiments provide an analysis which extends

15 the study of the system to include Lewis, X, Y, and precursor determinants in the entire human urinary tract.

20 One skilled in the art will recognize that the analysis of human urinary tract tissue is applicable to any tissue system for the determination of expression of blood group antigens. For example, the terms "secretor" and "non-secretor" are used to define individuals who do or do

25 not secrete A,B, or H antigens in saliva. "Secretors" produce Le<sup>b</sup> and Y antigens, whereas non-secretors produce Le<sup>a</sup> and X. Watkins, Science 152:172 (1966). As "secretor" or non-secretor" status, as well as changes in this status,

30 is considered indicative of cancer susceptibility or onset, the monoclonal antibody panels, and methods described herein, are useful in cancer diagnosis.

1

5 MATERIAL AND METHODSTissues

Human fetal tissues ranging from 12 to 14 weeks of gestational age were obtained from elective abortions.

10 Human normal adult tissues were obtained at autopsy within 9 hours post-mortem or from surgical pathology specimens within 1-2 hr of resection. Fresh tissues were fixed in 10% formaldehyde in phosphate buffered saline (PBS) (pH 7.4),

15 and embedded in paraffin. Alternatively, tissues were snap-frozen in isopentane precooled in liquid nitrogen, embedded in OCT compound in cryomolds and stored at -70°C until needed. Two fetal specimens containing kidney and

20 ureter were studied, one expressing A group and the other H group antigens. The adult kidney, ureter, and/or bladder tissues chosen for the present study included samples from 3 group O, 3 group A, 2 group B, and 2 group AB individuals.

25 The blood group of individuals from whom the specimens were derived was correlated with the immunohistological patterns of reactivities.

30 Reagents.

Purified agglutinin I from Ulex europaeus at 4 ug/ml served to identify the H-antigen. Mouse monoclonal antibody H 29-36 recognizes A antigen (all variants), mAb S8

35 detects B-antigen, mAbs T-174, T-218, P-12 and F-3 with

1

5 specificities for Le<sup>a</sup>, Le<sup>b</sup>, X and Y antigens, respectively,  
were also used. Finally, mAb K-21 detects precursor type I  
chain antigen. The antibodies were used as undiluted  
10 culture supernatants, or after purification from mouse  
ascities fluid (1:250 dilution).

Immunohistochemistry.

15 a) Indirect Immunofluorescence: Frozen tissues (4  
to 8 microns) were cut using a cryostat with a microtome.  
Cryostat-cut sections were used unfixed or fixed for 10  
minutes with either 1% formalin in PBS or cold acetone.  
Tissue sections were washed several times in PBS and rinsed  
20 in 2% bovine serum albumin in PBS (BSA-PBS). They were then  
incubated in a wet chamber with primary antibodies for 1  
hour at room temperature, the titration and appropriate  
dilution having been previously established. Sections were  
25 washed with PBS and incubated for 45 minutes with secondary  
fluoresceinated antibodies, which have also been previously  
titrated for optimal dilution (usually 1:40 in BSA-PBS).  
Tissue sections were washed extensively in PBS with the  
30 creation of turbulence, using a magnetic stirring plate, wet  
mounted in 90% glycerol in PBS, and examined with a  
fluorescence microscope equipped with epifluorescence, using  
a 100 watt mercury lamp.

35

1

5                   b) Immunoperoxidase: Formalin-fixed and  
paraffin-embedded sections were deparaffinized for this  
technique. Sections were treated for 30 minutes in 1%  
hydrogen peroxide in PBS in order to remove endogenous  
10 peroxidase activity (no staining was observed when 1%  
periodic acid was used instead of 1% hydrogen peroxide).  
Tissue sections were washed several times in PBS, and then  
incubated with the appropriate suppressor serum for 20  
15 minutes. Suppressor serum was drained off and sections were  
incubated with appropriately diluted primary antibody  
overnight at 4°C. Both peroxidase-antiperoxidase and  
avidin-biotin methods were used in these experiment. The  
20 secondary antibodies were horseradish peroxidase conjugated  
or biotinylated and they were incubated on sections for 1  
hr. Sections were then washed several times in PBS, and  
rinsed with 0.05M Tris buffer, 0.1 M NaCl, at pH 8. For the  
25 final reaction diaminobenzidine (DAB) was used as chromogen,  
and the peroxidase reaction was performed by incubating  
tissue sections for 6 to 12 minutes with 5 mg of DAB  
tetrahydrochloride in 100 ml of tris buffer containing 100  
30 ug of 0.3% hydrogen peroxide. Sections were washed with  
distilled water, counterstained with hematoxylin, and  
mounted with permount.

35

1  
5 Fresh frozen tissue sections were also used for this method. In this case, antibodies were incubated for 1 hr, and the other steps were similar to those described above for paraffin-embedded tissue sections.

10 c) Method for staining with lectin: The lectin Ulex europeus was incubated for 2 hours at room temperature, followed by rabbit anti-Ulex lectin antibody at a dilution  
15 of 1:1000 overnight at 4°C. Immunoperoxidase methods were performed as described above using biotinylated goat anti-rabbit immunoglobulins as secondary reagent.

20 d) Controls: Frozen and paraffin-embedded tissues expressing the appropriate blood group antigen served for titration of the reagents as well as positive and negative controls. Negative controls included substitution  
25 of the primary antibody by another antibody of the same species and isotype, or with PBS alone.

#### RESULTS

30 Table I summarizes the derivation of the panel of mouse mAb, their immunoglobulin subtype, and the characteristics of the blood group antigens detected. Table  
35 II summarizes the immunoreactivities of these antibodies on sections of normal adult kidney, ureter and urinary bladder.

1

5 Table III summarizes the immunoreactivities of this panel of antibodies on sections of fetal kidney and ureter. Figures 1 and 2 illustrate the immunohistological staining patterns of these mAb with normal human adult kidney and urothelium.

10

Blood group reactivities in adult tissues

Purified agglutinin I from Ulex europeus was used to identify the H-antigen. Expression of H-antigen was observed in endothelial cells and erythrocytes of all specimens studied. In the kidneys from A, B, or H individuals, H-antigen was found in the capillary network of glomeruli and also in the epithelial cells of collecting tubules with an homogeneous pattern of staining; the staining was weaker in the samples from AB individuals. All urothelial specimens expressed H-antigen throughout the mucosa, with an intense immunostaining of the basal layers. The rest of the nephron and connective tissue were negative for Ulex reactivity in all individuals tested.

20

Anti-A (H 29-36) and anti-B (S8) antibodies reacted only with the appropriately matched tissue specimens, that is from blood group A positive and B positive individuals, respectively. In each case, endothelial cells and erythrocytes were found to stain with the corresponding antibody. In the kidney, A and/or B antigens were found in glomerular and peritubular

30

35

1

5 capillaries and epithelial cells of collecting tubules (Fig. 1A). Urothelium was immunoreactive throughout, with some variation in staining intensity and usually greater reactivity in luminal cells (Fig. 1B).

10

Lewis antigens were expressed on the nephron with distinct patterns of reactivities in the adult kidney. Lewis<sup>a</sup> (T174) was generally observed in the adult kidney in the epithelial cells of collecting and distal tubules (Fig. 2A), and in one case (AB specimen) a faint staining of the proximal tubules and portions of the Henle's loop was observed. In urothelium, Lewis<sup>a</sup> was found to be consistently positive in the superficial epithelial cell layers (Fig. 2B) and weak or absent in deeper cell layers, though in one sample (H specimen) there was positive staining through all layers of the epithelium. Lewis<sup>b</sup> (T218) was found in rare collecting ducts and sometimes single cells of the collecting duct of adult kidneys, but with intense reactivity (Fig. 2C). In urothelium, Lewis<sup>b</sup> was expressed mainly in the basal and suprabasal cell layers (Fig. 2D), or in some specimens throughout the entire epithelium with increased intensity in the basal cells. X determinant (P12) was detected in polymorphonuclear leukocytes. In adult kidney, X antigen was positive in epithelial cells of proximal tubules, portions of the loop of Henle (Fig. 2E),

35

1

5 and in one case (A specimen) faint staining of distal and  
collecting tubular epithelial cells was observed. The  
reactivity of anti-X in urothelium was consistently intense  
in the umbrella cells, with only weak and variable staining  
10 of intermediate cell layers (Fig. 2F). Y determinant (F3)  
was detected in endothelial cells and erythrocytes. In  
adult kidney, endothelial cells of capillaries in the  
glomeruli were immunoreactive, as were the epithelial cells  
15 of collecting ducts (Fig. 2G). The reactivity with adult  
urothelium was intense and Y antigen was seen as an  
homogeneous pattern in the entire mucosa, with increased  
staining of basal and suprabasal cells (Fig. 2H).

20

Finally, the precursor type I chain (K21) was  
found in the adult kidney on occasional epithelial cells of  
the distal tubules and collecting ducts (Fig. 1C), and it  
25 was noted detected in any specimen of adult ureter or  
urinary bladder studied (Fig. 1D).

#### Blood group reactivities in fetal tissue

30

No differences were observed in the immunostaining  
patterns of fetal tissues when compared with those of the  
adult with the reagents detecting H, A, and B antigens.

35



1

5 In general, the expression of Le<sup>a</sup>, Le<sup>b</sup>, X, Y and  
precursor antigens in the fetal urinary tract resembled that  
of the adult although there were some significant  
differences (Table II). Lewis<sup>a</sup> staining was very strong in  
10 the collecting ducts of the fetal kidney, while staining of  
urothelium was weaker in fetus than in adult. X antigen  
reactivity in the kidney was similar to the adult; the fetal  
urothelium was more strongly reactive than adult. Le<sup>b</sup> and Y  
antigen reactivities were very similar in adult and fetal  
15 tissues throughout the urinary system.

Finally, in contrast to adult tissues, fetal  
kidney and ureter expressed precursor type I antigen as  
20 strong and homogeneous staining of epithelial cells of  
collecting tubules and urothelium.

One skilled in the art will see the applicability  
25 of the monoclonal antibodies and panels of monoclonal  
antibodies described herein. For example, tissues and  
organs may be typed according to expression of blood group  
antigens to determine if they may be used in  
30 transplantation. The blood group antigen profile of the  
organ and/or tissue to be transplanted is compared to the  
blood group antigen profile of the intended recipient to  
determine if cross reactivity is to be expected. Blood  
35 typing can be done also, again using the antibodies and

1

antibody panels described herein. A sample of blood is  
5 contacted to an antibody or panel of antibodies under  
conditions which favor complexing of the blood group  
antigens present with the known antibodies. By examining  
the pattern of complexing, one therefore determines the  
10 presence, or lack thereof, of particular blood group  
antigens. In typing organs and tissues, a similar practice  
is used (i.e., a cell sample or tissue sample is used for  
the contacting to the antibodies.

15

One skilled in the art will also see that in  
cancer patients, blood group antigen expression often  
changes, and different changes are characteristic of  
20 particular cancerous conditions. Thus, when cancer is  
suspected, a patient's blood, body secretions, or samples of  
tissue, are assayed using the monoclonal antibodies and  
antibody panels of this invention. Following the contacting  
25 method described supra, a pattern of blood group antigen  
expression is obtained, which is compared to an individual's  
normal blood group antigen expression panel, or "blood  
type." A diagnosis can then be made regarding the patient's  
30 condition with respect to cancer.

35

1

Hence, by using a panel of monoclonal antibodies  
5 including at least one antibody from the group consisting of  
H 29-36, S-8 T-174, T-218, P-12, F-3, and K-21, the  
aforementioned aspects of this invention are accomplished.  
As will be seen, supra, this invention is particularly  
10 useful in diagnosing cancer of the larynx, respiratory  
tract, and urinary bladder or urinary tract cancer.

One skilled in the art will see also that the  
15 antibodies and panel of antibodies described herein can be  
obtained in the form of kits, wherein different samples of  
monoclonal antibodies are separately packaged, such that  
individual antibodies, or the entire panel may be used, as  
20 desired. The individual samples allow one to perform  
sequential testing, for example.

While there have been described what are at  
25 present considered to be the preferred embodiments of this  
invention, it will be obvious to those skilled in the art  
that various changes and modifications may be made therein  
without departing from the invention, and it is, therefore,  
30 aimed to cover all such changes and modifications as fall  
within the true spirit and scope of the invention.

35

0232706

SK 375

TABLE II. IMMUNOREACTIVITY OF A PANEL OF ANTIBODIES DETECTING BLOOD GROUP ANTIGENS IN ADULT KIDNEY, URETER AND URINARY BLADDER.

BLOOD TYPE	TISSUES	ANTIGEN EXPRESSION							
		H	A	B	LeA	LeB	X	Y	PS
H(D)	GLOMERULUS	0	0	0	0	0	0	0	0
	PROXIMAL TUB.	0	0	0	0	0	0	0	0
	LOOP HENLE	0	0	0	0	0	0	0	0
	DISTAL TUB.	0	0	0	0	0	0	0	0
	COLLECTING TUB.	0	0	0	0	0	0	0	0
	UROTHELIUM	0	0	0	0/0	0/0	0	0	0
A	GLOMERULUS	0	0	0	0	0	0	0	0
	PROXIMAL TUB.	0	0	0	0	0	0	0	0
	LOOP HENLE	0	0	0	0	0	0	0	0
	DISTAL TUB.	0	0	0	0	0	0	0	0
	COLLECTING TUB.	0	0	0	0	0	0	0	0/0
	UROTHELIUM	0	0	0	0	0	0	0	0
B	GLOMERULUS	0	0	0	0	0	0	0	0
	PROXIMAL TUB.	0	0	0	0	0	0	0	0
	LOOP HENLE	0	0	0	0	0	0	0	0
	DISTAL TUB.	0	0	0	0	0	0	0	0
	COLLECTING TUB.	0	0	0	0	0	0	0	0
	UROTHELIUM	0	0	0	0	0	0	0	0
AB	GLOMERULUS	0	0	0	0	0	0	0	0
	PROXIMAL TUB.	0	0	0	0	0	0	0	0
	LOOP HENLE	0	0	0	0	0	0	0	0
	DISTAL TUB.	0	0	0	0	0	0	0	0
	COLLECTING TUB.	0	0	0	0	0	0	0	0
	UROTHELIUM	0	0	0	0	0	0	0	0

PS=PRECURSOR STRUCTURE ; IMMUNOREACTIVITIES: 0-HOMOGENEOUS STAINING,  
 0=OCCASIONAL CELLS, 0=LUMINAL SIDE POSITIVITY, 0=BASEL SIDE POSITIVITY.

TABLE III. IMMUNOREACTIVITIES OF A PANEL OF ANTIBODIES DETECTING BLOOD GROUP ANTIGENS IN FETAL KIDNEY AND URETER.

TISSUES	ANTIGEN EXPRESSION							
	H	A	B	Le <sup>A</sup>	Le <sup>B</sup>	X	Y	PS
GLOMERULUS	0	0	0	0	0	0	0	0
PROXIMAL TUB.	0	0	0	0	0	0	0	0
LOOP OF HENLE	0	0	0	0	0	0	0	0
DISTAL TUB.	0	0	0	0	0	0	0	0
COLLECTING TUB.	0	0	0	0	0	0	0	0
UROTHELIUM	0	0	0	0	0	0	0	0

PS = PRECURSOR STRUCTURE; IMMUNOREACTIVITIES: 0=HOMOGENEOUS STAINING,  
 0=OCCASIONAL CELLS, 0=LUMINAL SIDE POSITIVITY, 0=BASAL SIDE POSITIVITY.  
 TUB.=TUBULES.

1

5     What is Claimed:

1.     Panel for identifying blood group antigens  
          comprising at least one monoclonal antibody which  
10       specifically binds to a blood group antigen.
2.     Panel of Claim 1, comprising at least one monoclonal  
          antibody from the group consisting of H-29-36, S-8,  
15       T-174, T-218, P-12, F-3 and K-21.
3.     Panel of claim 1, comprising monoclonal antibodies H-  
          29-36, S-8, T-174, T-218, P-12, F-3, and K-21.  
20
4.     A method of typing blood comprising contacting a blood  
          sample to a panel of antibodies which specifically bind  
          to blood group antigens under conditions favoring  
25       formation of antibody-antigen complexes, and  
          determining which antibodies specifically bind to  
          antigens in said sample.
- 30     5.     A method as in Claim 4, wherein said method is used to  
          determine suitability of organ or tissue transplants.

35

- 1 6. A method of diagnosing disease associated with  
secretion of blood group antigens comprising contacting  
a tissue sample to a panel of monoclonal antibodies  
5 which specifically bind to different blood group  
antigens under conditions favoring formation of  
antibody-antigen complexes, determining presence of  
said complexes so as to determine presence of said  
10 antigens, and comparing antigen presence to antigen  
patterns associated with a particular disease.
7. Method of Claim 6, wherein said disease is cancer.  
15
8. Method of Claim 7, wherein said cancer is larynx,  
respiratory tract, urinary bladder or urinary tract  
cancer.  
20
9. A kit for use in determining the presence of blood  
group antigens comprising separate portion of at least  
one monoclonal antibody which specifically binds to a  
25 blood group antigen.
10. Kit of Claim 10, comprising at least one monoclonal  
antibody from the group consisting of H 29-36, S-8,  
30 T-174, T-218, P-12, F-3, and K-21.

FIGURE 1

