

This Page Is Inserted by IFW Operations
and is not a part of the Official Record

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images may include (but are not limited to):

- BLACK BORDERS
- TEXT CUT OFF AT TOP, BOTTOM OR SIDES
- FADED TEXT
- ILLEGIBLE TEXT
- SKEWED/SLANTED IMAGES
- COLORED PHOTOS
- BLACK OR VERY BLACK AND WHITE DARK PHOTOS
- GRAY SCALE DOCUMENTS

IMAGES ARE BEST AVAILABLE COPY.

**As rescanning documents *will not* correct images,
please do not report the images to the
Image Problem Mailbox.**

⑬



Europäisches Patentamt
European Patent Office
Office européen des brevets

⑪ Publication number:

0 232 706
A2

⑫

EUROPEAN PATENT APPLICATION

⑲ Application number: **87100099.8**

⑳ Date of filing: **07.01.87**

⑤① Int. Cl.: **G 01 N 33/80, G 01 N 33/52,**
G 01 N 33/574
// G01N33/577

③① Priority: **10.01.86 US 817874**

④③ Date of publication of application: **19.08.87**
Bulletin 87/34

⑥④ Designated Contracting States: **DE FR GB**

⑦① Applicant: **Sloan-Kettering Institute For Cancer Research, 1275 York Avenue, New York New York 10021 (US)**

⑦② Inventor: **Cordon-Cardo, Carlos, 430 East 67th Street, New York New York 10021 (US)**
Inventor: **Lloyd, Kenneth O., 4525 Henry Hudson Parkway West, Bronx New York 10471 (US)**
Inventor: **Finstad, Connie L., 504 East 81st Street, New York New York 10028 (US)**
Inventor: **Mc Groarty, M.E., 16 Franklin Place, Great Neck New York 10021 (US)**
Inventor: **Bander, Neil H., 1385 York Avenue, New York New York 10021 (US)**
Inventor: **Old, Lloyd J., 600 West End Avenue, New York New York 10024 (US)**
Inventor: **Melamed, Myron R., 35 Braeside Lane, Dobbs Ferry New York 10522 (US)**

⑦④ Representative: **Patentanwälte Schulze Horn und Hoffmeister, Goldstrasse 36, D-4400 Münster (DE)**

⑤④ **Blood group antigen panel.**

⑤⑦ Monoclonal antibodies specific for human blood group antigens are useful in cancer diagnosis. It has been found that cancer patients often either express antigens of a blood type different from their normal blood type or cease to produce normal blood group antigens in cancerous tissues. By examining tissue samples with the monoclonal antibodies described herein, it becomes possible to diagnose cancer.

EP 0 232 706 A2

1

FIELD OF THE INVENTION

5 This invention relates to the use of monoclonal
antibodies in determining the presence of particular
antigens. Applications of this invention include, but are
not limited to, diagnosis of disease, including cancer cell
10 typing or classification and identification of precancerous
lesions.

PRIOR ART

15 Blood group antigens are carbohydrate determinants
which are typically found on erythrocytes, certain
epithelial tissues, and in body secretions. They are formed
by the sequential addition of saccharides to carbohydrate
20 side chains of lipids and proteins. Hakomori, Seminars in
Hematology 18:39 (1981). Genes control synthesis of these
structures, as well as their expression in secretions, and
on cell types other than erythrocytes (red blood cells).

25 The "A", "B" and "H" blood group antigens are
known, at least indirectly, as identifying blood "type".
Presence of "H" antigen only is characteristic of "Type O"
30 blood, whereas presence of antigen "A" and "B" in the same
sample is characteristic of type "AB". Presence of "A" or
"B" antigen is characteristic of, respectively, type A or
type B blood. Lewis antigens, i.e., Le^a and Le^b, are
35 typically found in plasma, secretions, and secretory
epithelia.

1

These antigens are characteristic of other
5 conditions, and determination of their presence is useful in
areas other than blood typing. For example, Emmott, et al.,
J. Urol. 121:37 (1979) have found that the antigens of the
ABH system which are usually present in normal urinary
10 bladder tissue, are absent in urinary bladder tumors.
Additional studies have shown that loss of these antigens is
an early event in malignant transformation. Liss, et al.,
Am. J. Clin. Pathol. 68:372(1977); (larynx carcinoma);
15 Weinstein, et al., Cancer 43:661 (1979) (urinary bladder
carcinoma). In patients with epithelial cancers, especially
colon carcinomas, elevated levels of Le^a and Le^b antigens
have been found. Koprowski, et al., Science 212:53 (1981).
20 Additionally, the presence of normally incompatible blood
groups in the same patient has been described in some cancer
patients. Hatton, et al., Biochem. Biophys. Acta 666:361
(1981).

25

Determination of a change in the amount of blood
group antigen, a sudden appearance, or disappearance, is
indicative of a pathological state. Hatton, et al., supra,
30 for example, found A antigen in type O cancer patients.
Emmol, Lin, and Koprowski, supra have all shown that
disappearance or appearance of antigens is typical of
cancer. Hence, it is desirable to have a method for
35 determining the presence of blood group antigens.

1

5 Several monoclonal antibodies, which are known to
be specific for blood group antigens have been used to form
a "panel" for determining particular antigens. H29-36
monoclonal antibody, for example, determines the presence of
10 all varieties of A antigen. See, e.g., Sakamoto, et al.,
(unpublished manuscript), copending U.S. Patent Application
Serial No. 474,415, "Monoclonal antibody S8 is known to
detect B-antigen," Ueda, et al., PNAS 78:5122 (1981).
15 Additionally, monoclonal antibodies T-174, T-218, P-12, and
F-3 are specific for Le^a, Le^b, X, and Y antigens,
respectively. Antibody K-21 detects precursor type antigen.
Rettig, et al., Cancer Res. 45:815 (1985), Lloyd, et al.,
20 Immunogenetics 17:537 (1983).

 These monoclonal antibodies all of which are
described more fully, infra, are used as part of a panel to
25 determine blood group antigens.

 Further details on the panel, its uses, and
further embodiments, are presented in the description which
30 follows.

BRIEF DESCRIPTION OF THE FIGURES

 Figure 1 shows pictorially the structures and
35 origins of Le^a, Le^b, H-1, H-2, X, Y, A and B antigens.

1

5

Figures 2 and 3 illustrate immunohistological staining patterns of monoclonal antibodies described herein, when applied to normal human adult kidney and urothelium.

10

DETAILED DESCRIPTION OF PREFERRED EMBODIMENTS

Antibodies

15

The hybridoma cell lines which produce the monoclonal antibodies of this invention have been deposited with the American Type Culture Collection, 12301 Parklawn Drive, Rockville Maryland 20852 and bear the following accession numbers:

20

	<u>Hybridoma</u>	<u>ATCC #</u>
	H 29-36	HB 8248
	S 8	
25	T 174	HB 8242
	T 218	HB 8249
	P 12	HB 8551
	F 3	HB 8217
30	K 21	HB 8549

35

Information on derivation of these hybridomas may be found in copending U.S. Application Serial No. 474,415 (H 29-36, T-174, T-218), Serial No. 297,814 (S-8); Serial No.

1

5 604,080 (P-12) and K-21) and Serial No. 470,815 (F-3). In
addition, the hybridomas are described in Ueda, et al., PNAS
78:5122 (1981) (S8); Rettig, et al., Cancer Res. 45:815
10 (1985) (P-12, and K-21) and Lloyd, et al., Immunogenetics
17:537 (1983) (F-3). The disclosures of all of these are
incorporated by reference herein.

In summary, the hybridomas are prepared following
15 the Kohler-Millstein method well known to the art, using, as
immunizing cell lines, the materials set forth in Table 1.

20

25

30

35

0232706

1107-01-87

SK 375

Table 1. Derivation and specificity of mouse monoclonal antibodies identifying Blood Group Antigens.

Antibody (Ig subclass)	Immunizing Cell Line	Blood Group Specificity.
K21 (μ)	Tera-1 Teratocarcinoma	Precursor (Type 1 Chain)
T174 (γ 1)	SK-CO-10 Colon Cancer	Le ^a (Type 1 Chain)
T218 (μ)	SK-CO-10 Colon Cancer	Le ^b (Type 1 Chain)
P12 (μ)	Fresh Human Placenta	X (Type 2 Chain)
F3 (μ)	SK-1W-3 Lung Cancer	Y (Type 2 Chain)
T36 (γ 3)	HT29 Colon Cancer	A (Type 1 and 2 Chains)
S8 (μ)	SK-RC-7 Renal Cancer	B (Type 2 Chain)

Note: Purified agglutinin I from Ulex Europeans (Vector Laboratories, Burlingame, CA) served to identify the H-antigen.

1
5 Because expression of ABO(H) antigens on the
urogenital system has been discussed in the literature,
e.g., Coon, et al., Am. J. Clin. Path. 76:163 (1981);
Szulman, J. Exp. Med. 111:785 (1960), and recent studies
10 have shown localization of Le^a and Le^b antigens in normal
adult urothelium; Juhl, J. Histochem. Cytochem. 33:309
(1985), the urinary tract was used as an exemplary system.
The following experiments provide an analysis which extends
15 the study of the system to include Lewis, X, Y, and
precursor determinants in the entire human urinary tract.

20 One skilled in the art will recognize that the
analysis of human urinary tract tissue is applicable to any
tissue system for the determination of expression of blood
group antigens. For example, the terms "secretor" and
"non-secretor" are used to define individuals who do or do
25 not secrete A,B, or H antigens in saliva. "Secretors"
produce Le^b and Y antigens, whereas non-secretors produce
Le^a and X. Watkins, Science 152:172 (1966). As "secretor"
or non-secretor" status, as well as changes in this status,
30 is considered indicative of cancer susceptibility or onset,
the monoclonal antibody panels, and methods described
herein, are useful in cancer diagnosis.

35

1

5

MATERIAL AND METHODSTissues

Human fetal tissues ranging from 12 to 14 weeks of gestational age were obtained from elective abortions.

Human normal adult tissues were obtained at autopsy within 9 hours post-mortem or from surgical pathology specimens within 1-2 hr of resection. Fresh tissues were fixed in 10% formaldehyde in phosphate buffered saline (PBS) (pH 7.4), and embedded in paraffin. Alternatively, tissues were snap-frozen in isopentane precooled in liquid nitrogen, embedded in OCT compound in cryomolds and stored at -70°C until needed. Two fetal specimens containing kidney and ureter were studied, one expressing A group and the other H group antigens. The adult kidney, ureter, and/or bladder tissues chosen for the present study included samples from 3 group O, 3 group A, 2 group B, and 2 group AB individuals. The blood group of individuals from whom the specimens were derived was correlated with the immunohistological patterns of reactivities.

30

Reagents.

Purified agglutinin I from Ulex europaeus at 4 ug/ml served to identify the H-antigen. Mouse monoclonal antibody H 29-36 recognizes A antigen (all variants), mAb S8 detects B-antigen, mAbs T-174, T-218, P-12 and F-3 with

35

1

5 specificities for Le^a, Le^b, X and Y antigens, respectively,
were also used. Finally, mAb K-21 detects precursor type I
chain antigen. The antibodies were used as undiluted
10 culture supernatants, or after purification from mouse
ascities fluid (1:250 dilution).

Immunohistochemistry.

a) Indirect Immunofluorescence: Frozen tissues (4
15 to 8 microns) were cut using a cryostat with a microtome.
Cryostat-cut sections were used unfixed or fixed for 10
minutes with either 1% formalin in PBS or cold acetone.
Tissue sections were washed several times in PBS and rinsed
20 in 2% bovine serum albumin in PBS (BSA-PBS). They were then
incubated in a wet chamber with primary antibodies for 1
hour at room temperature, the titration and appropriate
dilution having been previously established. Sections were
25 washed with PBS and incubated for 45 minutes with secondary
fluoresceinated antibodies, which have also been previously
titrated for optimal dilution (usually 1:40 in BSA-PBS).
Tissue sections were washed extensively in PBS with the
30 creation of turbulence, using a magnetic stirring plate, wet
mounted in 90% glycerol in PBS, and examined with a
fluorescence microscope equipped with epifluorescence, using
a 100 watt mercury lamp.

35

1

5 b) Immunoperoxidase: Formalin-fixed and
paraffin-embedded sections were deparaffinized for this
technique. Sections were treated for 30 minutes in 1%
hydrogen peroxide in PBS in order to remove endogenous
10 peroxidase activity (no staining was observed when 1%
periodic acid was used instead of 1% hydrogen peroxide).
Tissue sections were washed several times in PBS, and then
incubated with the appropriate suppressor serum for 20
15 minutes. Suppressor serum was drained off and sections were
incubated with appropriately diluted primary antibody
overnight at 4°C. Both peroxidase-antiperoxidase and
avidin-biotin methods were used in these experiment. The
20 secondary antibodies were horseradish peroxidase conjugated
or biotinylated and they were incubated on sections for 1
hr. Sections were then washed several times in PBS, and
rinsed with 0.05M Tris buffer, 0.1 M NaCl, at pH 8. For the
25 final reaction diaminobenzidine (DAB) was used as chromogen,
and the peroxidase reaction was performed by incubating
tissue sections for 6 to 12 minutes with 5 mg of DAB
tetrahydrochloride in 100 ml of tris buffer containing 100
30 ug of 0.3% hydrogen peroxide. Sections were washed with
distilled water, counterstained with hematoxylin, and
mounted with permount.

35

1
5 Fresh frozen tissue sections were also used for this method. In this case, antibodies were incubated for 1 hr, and the other steps were similar to those described above for paraffin-embedded tissue sections.

10 c) Method for staining with lectin: The lectin Ulex europeus was incubated for 2 hours at room temperature, followed by rabbit anti-Ulex lectin antibody at a dilution
15 of 1:1000 overnight at 4°C. Immunoperoxidase methods were performed as described above using biotinylated goat anti-rabbit immunoglobulins as secondary reagent.

20 d) Controls: Frozen and paraffin-embedded tissues expressing the appropriate blood group antigen served for titration of the reagents as well as positive and negative controls. Negative controls included substitution
25 of the primary antibody by another antibody of the same species and isotype, or with PBS alone.

RESULTS

30 Table I summarizes the derivation of the panel of mouse mAb, their immunoglobulin subtype, and the characteristics of the blood group antigens detected. Table
35 II summarizes the immunoreactivities of these antibodies on sections of normal adult kidney, ureter and urinary bladder.

1

5 Table III summarizes the immunoreactivities of this panel of
antibodies on sections of fetal kidney and ureter. Figures
1 and 2 illustrate the immunohistological staining patterns
of these mAb with normal human adult kidney and urothelium.

10

Blood group reactivities in adult tissues

Purified agglutinin I from Ulex europeus was used
to identify the H-antigen. Expression of H-antigen was
15 observed in endothelial cells and erythrocytes of all
specimens studied. In the kidneys from A, B, or H
individuals, H-antigen was found in the capillary network of
glomeruli and also in the epithelial cells of collecting
20 tubules with an homogeneous pattern of staining; the
staining was weaker in the samples from AB individuals. All
urothelial specimens expressed H-antigen throughout the
mucosa, with an intense immunostaining of the basal layers.
25 The rest of the nephron and connective tissue were negative
for Ulex reactivity in all individuals tested.

Anti-A (H 29-36) and anti-B (S8) antibodies
30 reacted only with the appropriately matched tissue
specimens, that is from blood group A positive and B
positive individuals, respectively. In each case,
endothelial cells and erythrocytes were found to stain with
35 the corresponding antibody. In the kidney, A and/or B
antigens were found in glomerular and peritubular

1

5 capillaries and epithelial cells of collecting tubules (Fig. 1A). Urothelium was immunoreactive throughout, with some variation in staining intensity and usually greater reactivity in luminal cells (Fig. 1B).

10

Lewis antigens were expressed on the nephron with distinct patterns of reactivities in the adult kidney. Lewis^a (T174) was generally observed in the adult kidney in the epithelial cells of collecting and distal tubules (Fig. 2A), and in one case (AB specimen) a faint staining of the proximal tubules and portions of the Henle's loop was observed. In urothelium, Lewis^a was found to be consistently positive in the superficial epithelial cell layers (Fig. 2B) and weak or absent in deeper cell layers, though in one sample (H specimen) there was positive staining through all layers of the epithelium. Lewis^b (T218) was found in rare collecting ducts and sometimes single cells of the collecting duct of adult kidneys, but with intense reactivity (Fig. 2C). In urothelium, Lewis^b was expressed mainly in the basal and suprabasal cell layers (Fig. 2D), or in some specimens throughout the entire epithelium with increased intensity in the basal cells. X determinant (P12) was detected in polymorphonuclear leukocytes. In adult kidney, X antigen was positive in epithelial cells of proximal tubules, portions of the loop of Henle (Fig. 2E),

35

1

5 and in one case (A specimen) faint staining of distal and
collecting tubular epithelial cells was observed. The
reactivity of anti-X in urothelium was consistently intense
in the umbrella cells, with only weak and variable staining
10 of intermediate cell layers (Fig. 2F). Y determinant (F3)
was detected in endothelial cells and erythrocytes. In
adult kidney, endothelial cells of capillaries in the
glomeruli were immunoreactive, as were the epithelial cells
15 of collecting ducts (Fig. 2G). The reactivity with adult
urothelium was intense and Y antigen was seen as an
homogeneous pattern in the entire mucosa, with increased
staining of basal and suprabasal cells (Fig. 2H).

20

Finally, the precursor type I chain (K21) was
found in the adult kidney on occasional epithelial cells of
the distal tubules and collecting ducts (Fig. 1C), and it
25 was noted detected in any specimen of adult ureter or
urinary bladder studied (Fig. 1D).

Blood group reactivities in fetal tissue

30 No differences were observed in the immunostaining
patterns of fetal tissues when compared with those of the
adult with the reagents detecting H, A, and B antigens.

35

1

5

10

15

In general, the expression of Le^a, Le^b, X, Y and precursor antigens in the fetal urinary tract resembled that of the adult although there were some significant differences (Table II). Lewis^a staining was very strong in the collecting ducts of the fetal kidney, while staining of urothelium was weaker in fetus than in adult. X antigen reactivity in the kidney was similar to the adult; the fetal urothelium was more strongly reactive than adult. Le^b and Y antigen reactivities were very similar in adult and fetal tissues throughout the urinary system.

20

Finally, in contrast to adult tissues, fetal kidney and ureter expressed precursor type I antigen as strong and homogeneous staining of epithelial cells of collecting tubules and urothelium.

25

30

35

One skilled in the art will see the applicability of the monoclonal antibodies and panels of monoclonal antibodies described herein. For example, tissues and organs may be typed according to expression of blood group antigens to determine if they may be used in transplantation. The blood group antigen profile of the organ and/or tissue to be transplanted is compared to the blood group antigen profile of the intended recipient to determine if cross reactivity is to be expected. Blood typing can be done also, again using the antibodies and

1

antibody panels described herein. A sample of blood is
5 contacted to an antibody or panel of antibodies under
conditions which favor complexing of the blood group
antigens present with the known antibodies. By examining
the pattern of complexing, one therefore determines the
10 presence, or lack thereof, of particular blood group
antigens. In typing organs and tissues, a similar practice
is used (i.e., a cell sample or tissue sample is used for
the contacting to the antibodies.

15

One skilled in the art will also see that in
cancer patients, blood group antigen expression often
changes, and different changes are characteristic of
20 particular cancerous conditions. Thus, when cancer is
suspected, a patient's blood, body secretions, or samples of
tissue, are assayed using the monoclonal antibodies and
antibody panels of this invention. Following the contacting
25 method described supra, a pattern of blood group antigen
expression is obtained, which is compared to an individual's
normal blood group antigen expression panel, or "blood
type." A diagnosis can then be made regarding the patient's
30 condition with respect to cancer.

35

1

Hence, by using a panel of monoclonal antibodies
5 including at least one antibody from the group consisting of
H 29-36, S-8 T-174, T-218, P-12, F-3, and K-21, the
aforementioned aspects of this invention are accomplished.
As will be seen, supra, this invention is particularly
10 useful in diagnosing cancer of the larynx, respiratory
tract, and urinary bladder or urinary tract cancer.

One skilled in the art will see also that the
15 antibodies and panel of antibodies described herein can be
obtained in the form of kits, wherein different samples of
monoclonal antibodies are separately packaged, such that
individual antibodies, or the entire panel may be used, as
20 desired. The individual samples allow one to perform
sequential testing, for example.

While there have been described what are at
25 present considered to be the preferred embodiments of this
invention, it will be obvious to those skilled in the art
that various changes and modifications may be made therein
without departing from the invention, and it is, therefore,
30 aimed to cover all such changes and modifications as fall
within the true spirit and scope of the invention.

35

0232706

SK 375

TABLE II. IMMUNOREACTIVITY OF A PANEL OF ANTIBODIES DETECTING BLOOD GROUP ANTIGENS IN ADULT KIDNEY, URETER AND URINARY BLADDER.

BLOOD TYPE	TISSUES	ANTIGEN EXPRESSION							
		H	A	B	LeA	LeB	X	Y	PS
H(D)	GLOMERULUS	0	0	0	0	0	0	0	0
	PROXIMAL TUB.	0	0	0	0	0	0	0	0
	LOOP HENLE	0	0	0	0	0	0	0	0
	DISTAL TUB.	0	0	0	0	0	0	0	0
	COLLECTING TUB.	0	0	0	0	0	0	0	0
	UROTHELIUM	0	0	0	0/0	0/0	0	0	0
A	GLOMERULUS	0	0	0	0	0	0	0	0
	PROXIMAL TUB.	0	0	0	0	0	0	0	0
	LOOP HENLE	0	0	0	0	0	0	0	0
	DISTAL TUB.	0	0	0	0	0	0	0	0
	COLLECTING TUB.	0	0	0	0	0	0	0	0/0
	UROTHELIUM	0	0	0	0	0	0	0	0
B	GLOMERULUS	0	0	0	0	0	0	0	0
	PROXIMAL TUB.	0	0	0	0	0	0	0	0
	LOOP HENLE	0	0	0	0	0	0	0	0
	DISTAL TUB.	0	0	0	0	0	0	0	0
	COLLECTING TUB.	0	0	0	0	0	0	0	0
	UROTHELIUM	0	0	0	0	0	0	0	0
AB	GLOMERULUS	0	0	0	0	0	0	0	0
	PROXIMAL TUB.	0	0	0	0	0	0	0	0
	LOOP HENLE	0	0	0	0	0	0	0	0
	DISTAL TUB.	0	0	0	0	0	0	0	0
	COLLECTING TUB.	0	0	0	0	0	0	0	0
	UROTHELIUM	0	0	0	0	0	0	0	0

PS=PRECURSOR STRUCTURE ; IMMUNOREACTIVITIES: 0-HOMOGENEOUS STAINING,
 0-OCCASIONAL CELLS, 0=LUMINAL SIDE POSITIVITY, 0=BASEL SIDE POSITIVITY.

TABLE III. IMMUNOREACTIVITIES OF A PANEL OF ANTIBODIES DETECTING BLOOD GROUP ANTIGENS IN FETAL KIDNEY AND URETER.

TISSUES	ANTIGEN EXPRESSION							
	H	A	B	Le ^A	Le ^B	X	Y	PS
GLOMERULUS	0	0	0	0	0	0	0	0
PROXIMAL TUB.	0	0	0	0	0	0	0	0
LOOP OF HENLE	0	0	0	0	0	0	0	0
DISTAL TUB.	0	0	0	0	0	0	0	0
COLLECTING TUB.	0	0	0	0	0	0	0	0
UROTHELIUM	0	0	0	0	0	0	0	0

PS = PRECURSOR STRUCTURE; IMMUNOREACTIVITIES: 0=HOMOGENEOUS STAINING,
 0=OCCASIONAL CELLS, 0=LUMINAL SIDE POSITIVITY, 0=BASAL SIDE POSITIVITY.
 TUB.=TUBULES.

1

5 What is Claimed:

1. Panel for identifying blood group antigens
 comprising at least one monoclonal antibody which
10 specifically binds to a blood group antigen.
2. Panel of Claim 1, comprising at least one monoclonal
 antibody from the group consisting of H-29-36, S-8,
15 T-174, T-218, P-12, F-3 and K-21.
3. Panel of claim 1, comprising monoclonal antibodies H-
 29-36, S-8, T-174, T-218, P-12, F-3, and K-21.
20
4. A method of typing blood comprising contacting a blood
 sample to a panel of antibodies which specifically bind
 to blood group antigens under conditions favoring
25 formation of antibody-antigen complexes, and
 determining which antibodies specifically bind to
 antigens in said sample.
- 30 5. A method as in Claim 4, wherein said method is used to
 determine suitability of organ or tissue transplants.

35

- 1 6. A method of diagnosing disease associated with
secretion of blood group antigens comprising contacting
a tissue sample to a panel of monoclonal antibodies
5 which specifically bind to different blood group
antigens under conditions favoring formation of
antibody-antigen complexes, determining presence of
said complexes so as to determine presence of said
10 antigens, and comparing antigen presence to antigen
patterns associated with a particular disease.
7. Method of Claim 6, wherein said disease is cancer.
15
8. Method of Claim 7, wherein said cancer is larynx,
respiratory tract, urinary bladder or urinary tract
cancer.
20
9. A kit for use in determining the presence of blood
group antigens comprising separate portion of at least
one monoclonal antibody which specifically binds to a
25 blood group antigen.
10. Kit of Claim 10, comprising at least one monoclonal
antibody from the group consisting of H 29-36, S-8,
30 T-174, T-218, P-12, F-3, and K-21.
- 35

FIGURE 1

